



Laser Safety Eyewear

Laser Safety Eyewear

Our general purpose laser safety eyewear range represents excellent value for money despite making no compromises on quality or performance. The polycarbonate filters are robust enough to withstand all but the highest laser powers and are provided in a wide range of different wavelength formats. With over 125 laser filters and easy-to-wear frames that fit any face, there's a pair of laser safety eyewear for every user and every application.

The unprotected human eye is extremely sensitive to laser radiation and can be permanently damaged from direct or reflected beams. Protective eyewear in the form of goggles, glasses, or shields provides the principal means to ensure against ocular injury, and must be worn at all times during laser operation of Class IIIb (CDRH), 3B (IEC) or higher lasers.

Selecting the right type of eyewear is critical in reducing the amount of incident light to safe levels, while transmitting sufficient light for good vision. The following flowchart summarizes the decision-making plan for selecting the right protective eyewear for your laser and application.

A custom solution can also be developed quickly and at low cost. The pricing is filter dependent and independent of the frame style.

Global Laser supply laser safety eyewear for medical, dental, military, aerospace, scientific, communications and industrial applications.

CE-certified Laser Safety Eyewear meet US and international laser safety standards, including ANSI Z136.1, EN207/EN208/EN60825 and Z87.1 for impact.



Choosing Eye Protection

- 1 Facts - Specifications of the Laser and conditions of use.**
 - Consult your laser's manufacturer's userguide for eyewear requirements.
 - Calculate OD and power density requirements based on wavelength, power in watts (or for pulsed systems, using wavelength, power in joules, pulse length in seconds and pulse repetition rate in hertz), using laser safety software such as Lazan, LaserSafePC, Easy Haz or the LIA's Laser Hazard Evaluator Software.
 - Consider: Are there engineering controls limiting exposure to the beam? Is partial beam visibility required for alignment of visible beams? Is protection required for intra-beam exposure or is protection primarily for diffuse or scattered energy? For medical applications, are there different eyewear considerations for the Clinician and patient (in terms of VLT-visible light transmittance, full orbital coverage, weight of eyewear)? Will filter colour / colour rendition affect use? Are there multiple laser systems in the area, or is the eyewear designated for a single system?
- 2 Filter - OD, damage threshold and Visible Light Transmittance (VLT) requirements.**
 - Make sure the filter will reduce possible energy exposure to below the Maximum Permissible Exposure (MPE).
 - Check the Photopic Visible Light Transmittance (VLT) of the filter. VLT is the percentage of visible light transmitted through a filter, calculated against the spectral sensitivity of the eye to daylight. The higher the better. VLTs below 20% should be used in well-illuminated working environments.
- 3 Frame - Style and mode of wearing.**
 - Rule #1: if the glasses are uncomfortable, users will be tempted to not wear them.
 - Rule #2: Vanity rules, even in the lab. Users will wear what they like.
 - Many frames are designed to fit-over prescription glasses. Some are universal, fitting well for those who do and those who do not wear prescription glasses.
 - Ensure that the selected frame is face-forming, well-fitting with no gaps. Models with side-shields increase ambient light, cut down on obstructed viewing and decrease the non-beam hazard of walking into a door.
 - Polymer filters are available in the most variety of frames, often with the widest field of view and full angular coverage.
- 4 Fit - adjustability, comfort, vanity.**
 - Repeat of Rule #1: If the goggles don't fit, users won't wear them.
 - Repeat of Rule #2: Users won't wear what doesn't fit well or what they don't like.
- 5 Factors - additional considerations.**
 - Eye protection is only effective when worn. It's of no use if it's sitting on the shelf.
 - If eye protection is too heavy, poorly fitted, poorly designed or the VLT is too low, users will make the wrong choice: not to wear it.
 - Risk assessment must be part of the equation, use engineering controls to reduce the risk.
 - Filter technical data, including batch data, absorption characteristics, test reports, CE certificates and documentation of conformity should be available upon request or online.
 - Consider the source. You only get two eyes.

Fitover Style



Small universal fit, comfortable over prescription glasses or alone, offering a full field of view

31 140mmx132mmx50mm


CE certified



Medium or **Large** universal fit, comfortable over prescription glasses or alone. Modern styling and comfort -fit temples.

51
53 141mmx145mmx45mm
146mmx145mmx50mm

CE certified



Medium or **Large** universal fit, comfortable over prescription glasses or alone, offering a full field of view. Adjustable temples.

36
38 142mmx135mmx49mm
144mmx145mmx50mm

CE certified



Medium universal fit, comfortable over prescription glasses or alone. Modern styling and comfort -fit temples

33 145mmx138mmx54mm

CE certified



Extra Large universal fit, comfortable over prescription glasses or alone, offering a full field of view

39 148mmx140mmx54mm


CE certified



Medium universal fit, in Retro styling. Insert for RX or secondary filter

55 145mmx142mmx45mm


CE certified



Medium or **Large** universal fit, comfortable over prescription glasses or alone, offering a full field of view

700
900 145mmx145mmx54mm
150mmx148mmx56mm

CE certified



Medium universal fit, in comfortable over prescription glasses or alone. Adjustable temples; Flip down RX or secondary filter

40 143mmx130mmx48mm

CE certified

www.globalasertech.com

Switchboard: +44 (0)1495 212213 | UK Sales Office +44 (0)1495 322350

Wraparound Style



Universal soft comfort fit goggle with insert for RX or secondary filter. Anti-fog venting and adjustable strap

50

133mmx51mm

CE certified



Medium modern style wraparound with wide field of view, comfort-fit temples and a faceform fit.

52

130mmx148mmx45mm

CE certified



Medium wraparound style frame, with adjustable temple lengths and adjustable temples angles. Comfort nose pads.

32

142mmx130mmx50mm

CE certified



Medium modern style wraparound with comfort-fit temples and nose pads, rimless edges and faceform fit. Insert for RX or secondary filters.

34

135mmx142mmx48mm

CE certified



Medium universal soft comfort fit goggle with Anti-fog venting and adjustable strap. Comfortable over prescription glasses or alone.

60

150mmx85mm

CE certified



Medium modern style wraparound with wide field of view, adjustable temple angles for faceform fit.

35

130mmx130mmx47mm

CE certified



Pediatric or Petite wraparound with slender faceform fit

KM

115mmx130mmx36mm

CE certified



Medium modern style wraparound with comfort-fit temples and nose pads for a faceform fit.

56

137mmx145mmx43mm

CE certified

www.globalasertech.com

Switchboard: +44 (0)1495 212213 | UK Sales Office +44 (0)1495 322350

Caution

Never look directly into the path of a laser. The laser safety eyewear offers protection against incidental exposure to specified beam energy only. Consult a laser safety officer, appropriate laser safety standards and/or laser system documentation to ensure correct eyewear and safe laser use. The use of incorrect eyewear may lead to serious personal and/or blindness.

Please note: Global Laser reserve the rights to change descriptions and specifications without notice.



For further information about any of our products please contact your local distributor or you can contact Global Laser in the UK. Your Local Distributor Is:

T: +44 (0)1495 212213
F: +44 (0)1495 214004
E: sales@globalasertech.com
www.globalasertech.com

Global Laser Ltd
Unit9-10
Roseheyworth Business Park
Abertillery, Gwent NP13 1SP UK